

Klear Dual Enzyme Block Kit

For Blocking Endogenous Peroxidase and Alkaline Phosphatase

| |
|----------------|
| Storage: 2-8°C |
|----------------|

Catalog No: E36-18 18 mL
 E36-100 100 mL

Intended Use:

Klear Dual Enzyme Block Kit is used to block both endogenous alkaline phosphatase and peroxidase enzymatic activity which can be used in paraffin-embedded tissue sections. **Klear Dual Enzyme Block Kit** is provided with the ready-to-use solution in a dropper bottle.

Kit Components:

| Product | Catalog No. | Size | Number of Slide Tests |
|-------------------------------|-------------|-------|-----------------------|
| Klear Dual Enzyme Block (RTU) | E36-18 | 18mL | 180 slides |
| Klear Dual Enzyme Block (RTU) | E36-100 | 100mL | 1200 slides |

Recommended Protocol:

This blocking reagent should be finished prior to application of primary antibody.

After deparaffinization:

- 1) Apply 2 drops (50 µL/ drop) of **Klear Dual Enzyme Block** to cover the tissue and incubate for at least 10 minutes at room temperature. Additional blocking time may be required depending on tissue thickness and sample penetrability.
- 2) Rinse well with distilled water.

Storage:

Store at 2-8°C. Return reagent back to 2-8°C after each use.

Precautions:

Please wear gloves and take other necessary precautions.

Remarks:

For research use only.