

## Normal Donkey Serum Blocking Buffer

(Normal Donkey Serum)

Storage: 2-8°C

Catalog No: E27-18      18mL  
E27-100      100mL

  

### Intended Use:

Normal Donkey Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is Donkey.

### Kit Components:

Catalog No.	Size
E27-18	18mL (Ready to use)
E27-100	100mL (Ready to use)

Solution contains 0.05% NaN<sub>3</sub>.

### Recommended Protocol:

1. Apply 100uL or enough volume of **Normal Donkey Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Donkey Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

### Precaution:

Wear gloves and take other necessary laboratory safety procedure.

### Storage:

Store at 2-8°C.

### Remarks:

For research use only.