

Normal Equine (Horse) Serum Blocking Buffer

Normal Equine (Horse) Serum

Storage: 2-8°C

Catalog No: E28-18	18mL	<input type="checkbox"/>
E28-100	100mL	<input type="checkbox"/>

Intended Use:

Normal Equine (Horse) Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is Equine (Horse).

Kit Components:

Catalog No.	Size
E28-18	18mL (Ready to use)
E28-100	100mL (Ready to use)

Solution contains 0.05% NaN₃.

Recommended Protocol:

1. Apply 100uL or enough volume of **Normal Equine (Horse) Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Equine (Horse) Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

Remarks:

For research use only.