

Normal Rabbit Serum Blocking Buffer

Normal Rabbit Serum

Storage: 2-8°C

Catalog No: E33-18 18mL
E33-100 100mL

Intended Use:

Normal Rabbit Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is Rabbit.

Kit Components:

Catalog No.	Size
E33-18	18mL (Ready to use)
E33-100	100mL (Ready to use)

Solution contains 0.05% NaN₃.

Recommended Protocol:

1. Apply 100uL or enough volume of **Normal Rabbit Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Rabbit Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

Remarks:

For research use only.